




T1 ACCELEROMETER/GYRO ASSEMBLY MOUNT for BioRID Model ARA-620

INTRODUCTION

The Model ARA-620 mount package is designed to allow measurements of T1 rotation in a BioRID dummy. It enables the IES-3103 Gyro sensors to be mounted to the T1 load cell using the existing accelerometer mounting holes. It is mounted on the opposite side of the spine from the accelerometer cube.



4	2	9000675	SCREW, SHCS #2-56 X 3/8
3	1	ARA-621	PLATE, GYRO-SENSOR
2	2	9000139	SCREW, SHCS #6-32 X 3/4
1	1	IES3103	3-AXIS GYRO SENSOR(REF)
ITEM	QTY	NAME	DESCRIPTION

Denton ATD, Inc. 10317 U.S. Highway 250 North, Milan, Ohio 44846-9570
Tel (419) 625-5200 * Fax (419) 625-5335 * email: info@dentonatd.com * www.dentonatd.com

Creating the New Standard in ATDs

Revision: 04/05/07